ADD

anterior distraction device

ADD*plus*

anterior distraction device plus

OP-TECHNIK
SURGICAL TECHNIQUE

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

The small pictures at the bottom of the page show the step-by-step application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-colored.

Art. Nr./Art. No.

Art. Nr./Art. No.

The implant ADD (anterior distraction device) is a vertebral body replacement for the anterior cervical and thoracic spine. ADD can also be used for intersomatic fusion.

ADD is used for the reconstruction of substance defects, support of the anterior spine as well as restoration of the spinal profile. An additional anterior or posterior stabilization by means of a rod/screw system or a plate system is necessary.

ADD*plus* (anterior distraction device plus) is also a vertebral body replacement which can be used for the anterior cervical spine and for the cervico-thoracic transition. The ADD*plus* is the result of the consistent and continuous development of our two systems ADD and osmium. The decisive advantage is the direct combination of vertebral body replacement and anterior plating, two systems in one – cage and fixation system.

ADD*plus* effects the reconstruction, support and stabilization of the anterior spine. The fixation on the vertebral bodies is effected by means of screws from the osmium system, with a monocortical as well as a bicortical screw anchoring being possible. The additional plates prevent a migration as well as possible sinking in of the implant. The combined implant thus provides a gain in additional safety. In comparison to a ventral plate the surface is reduced so that a improved bone fusion can be achieved.

Both implants, ADD and ADD*plus*, are continuously expandable in situ, thus an optimum bridging of the defect is guaranteed.

The large support surfaces with spikes guarantee secure support and anchoring.

The consolidation of introduced bone material is supported. The easy handling and few operating steps provide for short operating times and low blood loss for the patient.

The present surgical technique describes the procedure necessary for the application of the ADD and ADD*plus*. Also, the implants and the system instruments required for implantation are presented.

The surgical technique is not sufficient as the sole basis for the successful application of the ADD and ADD*plus*. It is recommended to study and learn the operating technique with and from an experienced surgeon.

For the application of the osmium system, a relevant description is additionally provided in our osmium surgical technique.

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Indications and contraindications

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Indications

ADD

Instabilities of different genesis after complete or incomplete corpectomy resulting from vertebral body destruction due to tumors, fractures or infections.

ADDplus

Necessity for multisegmental spondylodesis in the upper spine (C3-Th3). Instabilities of various genesis due to complete or incomplete corpectomy as a result of vertebral body destruction due to e.g. tumor or fracture.

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Contraindications

ADD und ADDplus must not be used in the case of acute infections, and proven incompatibility with regard to the metals used. In the case of an unfavourable medical or psychological general condition of the patient, which could be worsened further by the operation, a careful evaluation of advantages and disadvantages must be carried out.

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
ADD

ADD implants (CS 2250-X) are available in three diameters (12, 14 and 16 mm). The expansion ranges extend from 10 to 65 mm. In part, the implants are available with an angle of 6°. By turning the expansion ring the implant is spread open and the expansion height becomes fixed. The expansion height of the implant may be adjusted continuously as per the individual case.

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
ADDplus

ADDplus implants (CS 2253-X) are available in one diameter (12 mm), four expansion ranges and four angulations:

Distraktionshöhe	Winkelung
13–18 mm	0°
17–26 mm	6°
25–41 mm	12°
40–65 mm	18°

expansion height	angulation
13–18 mm	0°
17–26 mm	6°
25–41 mm	12°
40–65 mm	18°

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Locking screw

The locking screw (CS 2259) is used for the fixation of the expanded ADD or ADDplus implant.

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osmium screw

The osmium screw (CS 1300-X) is available in different lengths (14, 16, 18 mm) with a diameter of 5 mm. The length is easily recognizable according to the color: 14 mm golden, 16 mm green and 18 mm blue. Expansion of the screw is achieved by inserting the bolt. Overexpansion of the screw is not possible.

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
Cancellous bone screw

The cancellous bone screw (CS 1301-X) is used for the bicortical fixation. The diameter is 4 mm. The lengths of the screws are 12, 14, 16, 18, 20, 22, 24 and 26 mm.

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
Insertor for large heights
(larger than 16 mm)

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

Depth adjuster and
locking nut
for CS 1324-1

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Tap for cancellous bone
screws, Ø 4.0 mm,
for CS 1324

▲ CS 1328

Schraubendreher,
SW 2,5 mm

Screwdriver, hex 2.5 mm

▲ UT 1068-23

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Depth gauge,
length 230 mm

ADD and ADDplus surgical technique

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Preparation

The implantation of an ADD or ADDplus implant is demonstrated on a cervical model. Depending on actual segment height, an anterior or lateral approach can be selected for the use of ADD. ADDplus is used via an anterior approach. The affected vertebral body is exposed. Next, the adjacent vertebral discs and the vertebral body are fully or partially resected.

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Selection of the implant


To determine the appropriate size of ADD respectively ADDplus, the gap to be bridged is measured using the measuring caliper (CS 5788).

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
Determining the size of the implant

Using the measuring templates, the size of the required implant is determined. For ADD, a separate measuring template (CS 2266-3) is available. For ADDplus, a measuring template integrated into the layer for ADDplus implants (CS 2258-7) is available. The white section corresponds to the non-expanded implant. The shaded section marks the possible expansion height.

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
Pre-expansion of the implant

The implant ADD or ADDplus can already be expanded by hand, so that the height is within the area of the bone defect to be bridged.

ADD and ADDplus surgical technique

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Fitting the locking screw

The locking screw (CS 2259) can be pre-fitted outside the situs into a bore below the expansion ring. In the case of ADDplus the locking screw must be used on the front of the implant using the screwdriver (CS 1328). If inclined ADD implants are used, ensure that the locking screw is located at the approach end. The locking screw is not tightened.

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Fitting the inserter

The inserter (CS 2255-1 or CS 2255-2) is inserted into the implant bores. By tightening the knurled screws at the end of the instrument the implant is fixed to the instrument.

Important: Ensure that the locking screw is located between the arms of the inserter and that the expansion ring remains freely accessible.

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Insertion of the implant

The implant ADD or ADDplus is brought into the situs. In the case of ADDplus ensure that the plate marked with "top" points to the cranial end. This allows the angled end to better fit the contour of the lower end plate of the top vertebral body. Please note that recontouring of the implant, especially in the region of the plates must be desisted. Bending can lead to breakage of the implant.

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
Expansion

The expansion instrument (CS 2252) is inserted into a bore of the expansion ring of ADD or ADDplus. The cage is expanded by rotation of the ring in direction of the arrow.

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
Locking the implant ADD

After checking the correct position and the height of the ADD the locking screw is tightened using the screwdriver (CS 2261).

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Further steps are only required for the ADDplus instrumentation.

The ADDplus is secured by four monocortical osmium screws or bicortical cancellous bone screws from the osmium system. For additional information, the osmium surgical technique can be used.

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
Assembling the depth stop

The depth adjuster (CS 1324-2) with the locking nut (CS 1324-3) is screwed into the osmium depth stop (CS 1324-1). The drilling and tapping depth is determined by the x-ray images or CT/MRI images and adjusted on the depth stop. The depth stop is also used as a tissue protector. The drilling depth corresponds to the screw length.

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
Bore for osmium screw

The drill hole is prepared using the drill bit, Ø 2.5 mm (CS 1323) through the osmium depth stop (CS 1324). Drill under C-arm control, if necessary.

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
Tapping for osmium screw

The tap for the screw is prepared using the tap Ø 5 mm (CS 1325) through the osmium depth stop.

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

Preparing screw insertion

First remove the bolt of the osmium screw with the screwdriver (CS 1328). The attachment (CS 1310-2) is placed in the screwdriver for screw (CS 1310-1). By turning the attachment the osmium screw is locked. The prongs of the screwdriver must be fitted to the grooves of the screws.

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Implanting the osmium screws

Using the screwdriver (CS 1310-1) or the tightener for osmium screws (CS 1322) the screws can be further tightened before being expanded, in order to achieve an even support of the plates.

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Fixing the osmium screws

The screwdriver for osmium screw (CS 1310-1) without the attachment (CS 1310-2) is placed onto the screw and the bolt is inserted into the shaft of the instrument.

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

Fixing the osmium screws

The bolt is locked using the screwdriver, hex 2.5 mm (CS 1328) and the osmium screw is expanded.

Important: The instrument must remain engaged on the grooves of the osmium screw, in order to prevent that the screw corotates.

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Locking the implant ADDplus

After checking the correct positioning and height of the ADDplus the locking screw is tightened by using the screwdriver (CS 2261).

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Bicortical instrumentation with cancellous bone screws

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Preparing screw implantation

For the cancellous bone screws the same drill bit (CS 1323) and the depth stop (CS 1324) may be used. Determining the drill depth by using the depth stop is possible as well. The drilling depth corresponds to the screw length.

Important: Drill under C-arm control!

Tapping for cancellous bone screws

Use the same technique as with the monocortical screws. However, use the tap for cancellous bone screws, Ø 4 mm (CS 1326).

Important: It is optional to use the tap (CS 1326) before applying the screws. Use C-arm control!

Implanting the cancellous bone screws

The appropriate length of the cancellous bone screw may also be determined using the depth gauge (UT 1068-23). Inserting the cancellous bone screws by means of the screwdriver, hex 2.5 mm (CS 1328) under C-arm control.

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




Instruments	Product number
Screw holding forceps	CS 1219
Screwdriver for osmium screw	CS 1310-1
Attachment for screwdriver	CS 1310-2
osmium screw tightener	CS 1322
Drill bit , Ø 2.5 mm, for CS 1324	CS 1323
Depth stop	CS 1324-1
Depth adjuster for CS 1324-1	CS 1324-2
Locking nut for CS 1324-1	CS 1324-3
Tap for osmium screws , Ø 5 mm, for CS 1324	CS 1325
Tap for cancellous bone screws , Ø 4 mm, for CS 1324	CS 1326
Screwdriver , hex 2.5 mm	CS 1328
Expansion instrument ADD	CS 2252
Insertor for small heights	CS 2255-1
Insertor for large heights (larger than 16 mm)	CS 2255-2
Screwdriver , hex 2.5 mm	CS 2261
Measuring template for ADD	CS 2266-3
Impactor , height 8 mm, straight	CS 2274-2
Measuring caliper	CS 5788
Depth gauge , length 230 mm	UT 1068-23

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Further information are available here:

- www.ulrichmedical.com
- personal, per phone, fax, e-mail
- or at the next congress