# MENDEC<sup>®</sup> SpINE

High Barium Content Acrylic Resin

*for OUS market*

## TECRES<sup>®</sup>

*Advancing High Technology*

TECRES S.p.A.

Via A. Doria, 6 - 37066 Sommacampagna • Verona - ITALY

Telefono +39.045.9217311 - Fax +39.045.9217330

info@tecres.it • www.tecres.it

for Tecres  
internal use

CE  
0086

# MENDEC® Spine

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## PRESENTATION

MENDEC SPINE is an acrylic resin radiopaque for the filling of pathologic vertebral bodies. MENDEC SPINE is a high fluidity resin for syringe application.

MENDEC SPINE is a sterile and single-use medical device. The package consists of a container of powder and a vial of liquid, both sterilised using ethylene oxide. The liquid is sterilised through filtration, and the powder is sterilised with ethylene oxide.

Composition of MENDEC SPINE:

POWDER	20 g container of sterile powder
Polymethylmethacrylate	67,5% w/w
Barium sulphate	30,0% w/w
Benzoyl peroxide	2,5% w/w
LIQUID	9,4 g phial of sterile liquid
Methylmethacrylate	99,1% w/w
N,N-dimethyl-p-toluidine	0,9% w/w
Hydroquinone	75 ppm

## INDICATIONS

Filling of pathologic vertebral body in relation to osteoporotic compression fractures, metastasis, myeloma.

The use in children is advised only when it is believed impossible to operate through other forms of intervention.

## CONTRAINDICATIONS

Active infections

Emorrhagic diathesis

Extended vertebral destruction, with collapse superior to 2/3 of the standard thickness, destruction of the posterior wall with epidural extension of the pathologic tissue and clinical signs of medullar compression.

The use of MENDEC SPINE must be carefully considered in presence of hypersensitivity to monomer or to any of the other components of the resin.

## USE DURING PREGNANCY AND BREAST-FEEDING

There are no tests which demonstrate the utilisation safety of MENDEC SPINE during pregnancy or breast-feeding. MENDEC SPINE should not be used in the first three months of pregnancy; for the remaining gestation period, should only be used in life-endangering situations.

## PRECAUTIONS FOR USE

- Read this instruction booklet very carefully.
  - Make sure that the inner package is undamaged and that the components are integral, i.e. the powder should not present yellow or brown discolouring and the liquid should not be syrupy.
  - Temperature has a major effect on the preparation characteristics of any acrylic resin. Temperatures of more than 23°C for the product, the preparation accessories and the environment accelerate the various stages in the preparation procedure. Lower temperatures retard the preparation stages. **Before using MENDEC SPINE it is strongly advised to make sure that the package is stored at a temperature of 23°C ± 1°C for the previous 24 hours.**
  - Make sure that the preparation accessories are specifically compatible with the product.
  - Use syringes and needles of proven chemical compatibility with the resin components.
  - Do not open the phial of liquid over the mixing bowl to avoid the risk of glass fragments entering the dough.
  - Do not mix the resin in a flow of air since this will promote rapid evaporation of the liquid component and consequent variations in the performance of the resin.
  - NEVER add other substances or foreign bodies to the acrylic resin.
- To ensure safe application of MENDEC SPINE, surgeons must be aware of the properties of the product, its preparation and handling characteristics, applicational limitations and correct application technique. It is strongly recommended that the surgical team carry out practical trials prior to use in patients under the same instrumental and environmental conditions.

Clinical data demonstrate the need to maintain strictly aseptic surgical techniques. It is important to keep in mind that the possible deep infection of a surgical wound is a serious risk for the successful outcome of the surgical operation.

Prolonged respiration of the vapours may cause drowsiness: excessive exposure to monomer vapours may irritate the respiratory apparatus and the eyes. Avoid monomer contact with the skin and mucous membranes. Cases of contact dermatitis have been observed in susceptible subjects. It is therefore advisable to wear a second pair of surgical gloves and scrupulously observe the instructions for mixing the components in order to reduce the possibility of reactions caused by hypersensitivity. Do not allow the liquid component to come into direct contact with accessories in rubber or elastomers. The vapours of the liquid component may affect soft contact lenses.

The liquid component is flammable and volatile and for this reason the operating theatre must be correctly ventilated. The liquid component and/or its vapours must never be directly exposed to naked flames or heated materials.

**Caution:** Acrylic resins reach temperatures higher than physiological temperatures during the exothermic polymerisation reaction.

Following surgery, if any form of infection should arise, patients must immediately consult their doctors to reduce the risk of infection.

**PREPARATION AND APPLICATION INSTRUCTIONS**

Set the container of powder on a flat surface and open the cap. Open the vial and pour all of the liquid into the powder. Close the container and manually shake it until a liquid and homogenous mass is achieved. Set the container down again, remove the cap, draw the resin into a syringe and dispense it. If you are using a resin supply device with a luer-lok connection: open the resin container's luer top and connect it with the luer-lok connection in order to suck the resin from the resin supply device.

**Caution: never arbitrarily modify the ratios between the liquid and solid components.**

Do not re-sterilize and/or re-use the device. It is designed for single-use on a single patient.

Avoid dividing the product into two or more portions, in order to use it other times. This would be a re-utilization and may lead to an error in the proportion of the powder and liquid components. It can also cause a sterility loss.

Re-sterilization should not be carried out since it can cause infection risks for the patient. Re-sterilization can also alter the device's morphology and the mechanical features, causing a malfunction of the same with serious risks for the patient's health.

All residues must be considered waste surgical material and therefore should be eliminated at the end of the operation.

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The injection must be performed under continuous radioscopic control and must be done slowly, verifying the exit of the resin from the needle, its distribution inside the vertebral body, and checking for any possible extra-vertebral diffusion. If resin escapes from the vertebra, immediately discontinue the injection.

**TIMING FOR THE PREPARATION AND APPLICATION**

Times indicated in the instructions were obtained in laboratory by performing tests under temperature controlled conditions of  $23^{\circ}\text{C} \pm 1^{\circ}\text{C}$  and by using 20 ml plastic syringe.

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Operation:	Definition	Phase duration
Mixing	Mixing of the components	60 sec
Delivery device filling	The dough is transferred into the delivery device	60 sec
Waiting	The cement cannot be used	300 sec
Working	The cement can be delivered	600 sec
Hardening	The cement hardens increasing its viscosity and cannot be delivered anymore. Then the exothermic reaction takes place	360 sec

**If you are using an injection system that needs to be connected to an extension tube, let a few drops of resin pour out before connecting. If the resin is injected early it will result very liquid. If the resin is injected late it will result highly viscous.**

**WARNINGS**

- Store at a temperature below  $25^{\circ}\text{C}$  and away from light.
- Sterility is assured only if the unit container is not damaged or opened.
- Do not re-sterilise any of the components.
- Do not use the product if the powder has a yellowish or brownish colour or if the liquid is syrupy. These two conditions indicate that the product has not been stored correctly.

**NEGATIVE EFFECTS**

The following negative reactions may arise with the use of the acrylic resin, but they are not directly due to the acrylic resin itself.

The surgeon must be aware of these possible events and be ready to treat them if they appear.

- cardiac arrest
- myocardial infarct
- pulmonary embolism
- temporary drop in blood pressure
- haemorrhage and haematoma
- infection of surface/deep surgical lesion
- short-term irregularity of the cardiac rhythm
- sudden death
- diffusion of the resin outside the vertebral body: in the perivertebral veins (pulmonary embolism), in the epidural plex (mielopathy, radiculopathy), in the intervertebral disc.

**IMPORTANT INFORMATION FOR THE PHYSICIAN**

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**Shelf-life**

Shelf life of Mendec Spine is three years

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STERILE USING ASEPTIC PROCESSING TECHNIQUES

FLAMMABLE LIQUID

CONTENTS STERILE UNLESS PACKAGE IS DAMAGED OR OPENED

CAUTION: EXTERNAL SURFACE OF FIRST CONTAINER INSIDE IS NOT STERILE

USE BY

BATCH NUMBER

DO NOT RE-USE

DO NOT RESTERILISE

KEEP AWAY FROM SUNLIGHT

CONSULT INSTRUCTIONS FOR USE

# MENDEC® Spine

INSTRUCTIONS FOR USE  
To the Medical Staff's attention

## CAUTION

Federal law restricts this device to sale by or on the order of a physician. Users should be trained in use of percutaneous vertebral body injection procedures using bone cement for the treatment of vertebral pathologies and in the management of complications resulting from procedures.

Read this document carefully. Failure to follow instructions may cause malfunction and/or cause patient injury.

## DEVICE DESCRIPTION

MENDEC SPINE is a highly radiopaque acrylic resin of medium viscosity.

MENDEC SPINE is a sterile and single-use medical device. The package consists of a container of powder and a vial of liquid, both sterilized using ethylene oxide. The liquid is sterilised by means of filtration and the powder by means of ethylene oxide treatment.

Composition of MENDEC SPINE:

POWDER	20 g container of sterile powder
Polymethylmethacrylate	67,50% w/w
Barium sulphate	30,00% w/w
Benzoyl peroxide	2,50% w/w
LIQUID	9,4 g phial of sterile liquid
Methylmethacrylate	99,10% w/w
N,N-dimethyl-p-toluidine	0,90% w/w
Hydroquinone	75 ppm

## INDICATIONS

Mendec Spine is indicated for the treatment of pathological fractures of the vertebrae using a vertebroplasty or kyphoplasty procedure. Painful vertebral compression fractures may result from osteoporosis, benign lesions (hemangioma), and malignant lesions (metastatic cancers, myeloma).

## CONTRAINDICATIONS

- Active or incompletely treated infections at the site where the resin is to be applied
- Hemorrhagic diathesis
- Extended vertebral destruction, with collapse superior to 2/3 of the standard thickness, destruction of the posterior wall with epidural extension of the pathologic tissue and clinical signs of medullar compression
- Coagulation disorders or severe cardiopulmonary disease
- Spinal stenosis (>20% by retropulsed fragments)
- Compromise of the vertebral fractures due to posterior involvement
- Patient clearly improving on medical therapy
- Prophylaxis in metastatic or osteoporotic patients with no evidence of acute fracture
- Non-pathological, acute traumatic fractures of the vertebra
- Vertebra plana (collapse>90%)
- Compromise of the vertebral body or the walls of the pedicles
- Unstable vertebral fractures due to posterior involvement
- In patients previously sensitized to the monomer or any of the other components of the resin

## USE DURING PREGNANCY, BREAST-FEEDING AND CHILDREN

The safety and effectiveness of the Mendec Spine has not been studied in pregnant women, breast-feeding women or children.

## WARNINGS

- Read and understand instructions before using this product
- This product is intended for use only by clinicians with adequate training and experience in the percutaneous vertebral body injection procedures. The physician is advised to consult the medical literature regarding techniques, complications and hazards associated with the intended procedures.
- Tecres as a manufacturer does not recommend a surgical technique. It is the responsibility of the physician to determine the appropriateness of this product and specific technique for each patient
- Store at a temperature below 25°C and away from light.
- Sterile only if the package is unopened and undamaged.
- Do not use the product if the powder has a yellowish or brownish colour or if the liquid is syrupy. These two conditions indicate that the product has not been stored correctly.
- The users must have specific training and experience to be thoroughly familiar with the properties, handling characteristics and application of the product and percutaneous cement delivery. It is strongly recommended that the surgical team carry out practical trials prior to use in patients under the same environmental conditions using the same instruments for clinical use.

- Adverse patient reactions affecting the cardiovascular system have been associated with acrylic resins (bone cements): hypotensive reactions have occurred between 10 and 165 seconds. They have lasted from 30 seconds to 5 minutes. Some patients have progressed to cardiac arrest. Patients should be monitored carefully for any change in blood pressure during and immediately following the application of bone cement.
- The liquid component is flammable and volatile: the operating theatre must be correctly ventilated. The liquid component and/or its vapours must never be directly exposed to naked flames or heated materials.
- MMA can cause hypersensitivity reactions in susceptible persons, which may result in an anaphylactic response.
- The use of this product is not recommended in patients that do not exhibit a pathological condition, such as primary or secondary osteoporosis or a tumor, which would impair the ability of the patient to heal using conservative treatment methods.
- Adhere strictly to good surgical principles and techniques. Deep wound infection is a serious postoperative complication and may require total removal of the embedded cement. Deep wound infection may be latent and not manifest itself even for several years postoperatively.
- Follow handling, mixing and preparation instructions carefully.
- Do not re-sterilize and/or re-use the device. It is designed for single-use on a single patient. Re-sterilization should not be carried out since it can cause infection risks for the patient. Re-sterilization can also alter the morphology and the mechanical features of the device, causing a malfunction of the same with serious risks for the patient's health.
- Caution should be used to prevent excessive exposure to the concentrated vapors of the monomer, which may produce irritation of the respiratory tract, eyes and possibly liver damage.
- Do not use monomer if the liquid component shows any sign of thickening or premature polymerization.
- Do not allow the liquid component to contact rubber or latex gloves. The liquid component is a powerful lipid solvent. Should contact occur, the gloves may dissolve and tissue damage may occur. Wearing a second pair of gloves may diminish the possibility of hypersensitivity reactions.
- Do not allow personnel wearing contact lenses to be near or involved in mixing of the bone cement.
- Use appropriate imaging techniques to verify correct needle placement, absence of damage to surrounding structures, and appropriate location of the injected bone cement. Use imaging, such as fluoroscopy, to assess the ability of the vertebra to contain the injected bone cement.
- Cement leakage may cause tissue damage, nerve or circulatory problems, and other serious adverse events.
- Avoid over pressurization of the bone cement because this may lead to leakage of the bone cement beyond the site of its intended application. Cement leakage may cause tissue damage, and nerve or circulatory problems.
- Leaks can also occur when injecting if the needle is in a vein or if unseen microfractures are prevalent.
- If bone cement is seen outside of the vertebral body or in the circulatory system during the procedure, immediately stop the injection.
- Inadequate fixation or unanticipated postoperative events may affect the cement-bone interface and lead to micromotion of the cement against bone surface. A fibrous tissue layer may develop between the cement and the bone. Long term follow-up is advised for all patients on a regular scheduled basis.
- The completion of polymerization occurs in the patient and is an exothermic reaction with considerable liberation of heat. According to ISO 5833 standard, the temperature can be as high as  $90^{\circ}\text{C} \pm 5^{\circ}\text{C}$ . The heat produced in situ has a long term effect, and positioning of the patient should be maintained securely to obtain proper fixation. For proper fixation, 1-2 hours or longer may be required, as determined by the patient's medical condition and the attending physician.
- The long term effects of the bone cement in the spine have not yet been established.

#### PRECAUTIONS

- The product should not be used after the expiration date
- Avoid monomer contact with the skin and mucous membranes. Cases of contact dermatitis have been observed in susceptible subjects. It is therefore advisable to wear a second pair of surgical gloves and scrupulously observe the instructions for mixing the components in order to reduce the possibility of reactions caused by hypersensitivity.
- Temperature has a major effect on the preparation characteristics of any acrylic resin.
- Temperatures of more than  $23^{\circ}\text{C}$  for the product, the preparation accessories and the environment accelerate the various stages in the preparation procedure. Lower temperatures retard the preparation stages. **Before using MENDEC SPINE it is strongly advised to make sure that the package has been stored at a temperature of  $23^{\circ}\text{C} \pm 1^{\circ}\text{C}$  for the previous 24 hours.**
- Acrylic resins reach temperatures higher than physiological temperatures during the exothermic polymerization reaction. The released heat may damage bone or tissue surrounding the implant.
- Extrusion of the resin beyond the region of its intended application may occur.
- The polymer component should be disposed of in an authorized waste facility. The liquid component can be evaporated under a ventilated hood or absorbed by an inert material and transferred in a suitable container for disposal.
- Following surgery, patients should be instructed to contact their physicians if any signs or symptoms of infection occur.
- Adequately ventilate the operating room to eliminate as much monomer vapor as possible. The liquid monomer is highly volatile and flammable. Ignition of monomer fumes caused by the use of electrocautery devices in surgical sites near freshly implanted bone cements has been reported.

#### ADVERSE EVENTS

Serious adverse events, some with fatal outcome, associated with the use of acrylic bone cements for vertebroplasty include:

- cardiac arrest
- cardiac embolism
- cement extrusion into soft tissue
- cerebrovascular accident
- collapse of a vertebra adjacent to the one injected due to an osteoporotic disease

- compression of the spinal cord with paralysis or loss of feeling
- diffusion of the resin outside the vertebral body, such as the perivertebral veins (pulmonary embolism), the epidural plex (myelopathy, radiculopathy) or the intervertebral disc
- fracture of a pedicle
- haemorrhage and haematoma
- infection of surface/deep surgical lesion
- intercostal neuralgia
- myocardial infarction
- pneumonia
- pneumothorax
- pulmonary embolism
- rib fracture in patients with diffuse osteopenia, especially during thoracic vertebroplasty procedures, due to the significant downward force exerted during needle insertion
- short-term irregularity of the cardiac rhythm
- sudden death
- temporary drop in blood pressure

Although the majority of these adverse events present early within the post-operative period, there have been some reports of diagnoses beyond a year or more after the procedure.

Other reported adverse events for acrylic bone cements intended for vertebroplasty or kyphoplasty include:

Leakage of the bone cement beyond the site of its intended application with introduction into the vascular system resulting in embolism of the lung and/or heart or other clinical sequelae.

### **IMPORTANT PHYSICIAN INFORMATION**

- Percutaneous vertebroplasty procedures should only be performed in medical settings in which emergency decompressive surgery is available.
- The physician is responsible for any complication or harmful consequences which may result from an erroneous indication or operating technique, from inappropriate use of the material, or non-observation of the safety instructions that appear in the directions for use.
- Adverse reactions affecting the cardiovascular system have been reported following the application of polymethylmethacrylate bone cements. Recent data indicate that the monomer undergoes rapid hydrolysis to methacrylic acid, and that a significant fraction of the circulating methacrylic acid is in the form of circulating concentrations of methyl methacrylate/methacrylic acid although changes in blood pressure in response to these circulating concentrations have not been established.
- NEVER add other substances (including antibiotics) or foreign bodies to the acrylic resin. Moreover, never arbitrarily modify the ratios between the liquid and solid components.

The use of bone cement requires extremely close cooperation between the surgeon and the anaesthetist.

The anaesthetist must be told when the bone cement is implanted during the operation.

In some cases, events defined as bone cement implantation syndrome (BCIS) may occur, which are characterised by a series of clinical characteristics that include hypoxia, hypotension, cardiac arrhythmia, increased pulmonary vascular resistance and cardiac arrest, which must be controlled with the methods in use in modern anaesthesiology.

Even though these phenomena are rare, they are described in literature: percutaneous injection of bone cement can increase the intramedullary pressure and hence force both the bone marrow and the bone cement into the perivertebral venous system and these may migrate to the right side of the heart and into pulmonary circulation (pulmonary embolism) resulting in BCIS. BCIS has a very wide spectrum of severity and in many cases patients with pulmonary embolism detected during a chest X-ray remain asymptomatic.

### **INSTRUCTIONS FOR USE**

#### **Preparation**

Any **increase in temperature** of the working environment, the product and the preparation accessories over 23°C **reduces the waiting and hardening times** of the acrylic resin. Equally, **lower temperatures increase such times**.

Make sure that the inner package is undamaged and that the components are integral, i.e. the powder should not present yellow or brown discolouring and the liquid should not be syrupy.

Make sure that the preparation accessories are specifically compatible with the product.

Use accessories of proven chemical compatibility with the resin components.

#### **Mixing**

Set the container of the powder on a flat surface and open the cap.

Break open the vial and pour the liquid into the powder container. Do not open the vial of liquid over the powder container to avoid the risk of glass fragments entering the dough.

Close the container and manually shake it for 1 minute until a liquid and homogeneous mass is achieved.

**Avoid dividing the product into two or more portions, in order to use it other times. This would be a re-utilization and may lead to an error in the proportion of the powder and liquid components. It can also cause a sterility loss.**

Set the container down again, remove the cap and transfer the resin into the delivery apparatus. If you are using a resin supply device with a luer-lok connection: open the resin container's luer top and connect it with the luer-lok connection in order to suck the resin from the resin supply device. All residues must be considered waste surgical material and therefore should be eliminated at the end of the operation.

#### **Application**

Once the acrylic resin has reached the desired viscosity (waiting period) according to the delivery apparatus used, the injection may start: the injection must be performed under continuous radioscopic control and must be done slowly, verifying the exit of the resin inside the vertebral body and its distribution, and checking for any possible

extra-vertebral diffusion. If resin escapes from the vertebra, immediately discontinue the injection.

It is not advisable to use the product with needles with a caliber thinner than 15G.

Caution: never deliver acrylic resin into the vertebral body without the use of fluoroscopic guidance.

#### **TIMING FOR THE PREPARATION AND APPLICATION**

Times indicated in the instructions were obtained in the laboratory by performing tests under temperature controlled conditions of  $23^{\circ}\text{C} \pm 1$  ( $73^{\circ}\text{F} \pm 2$ ) and by using 20-ml plastic syringe.


Any **increase in temperature** of the working environment, the product and the preparation accessories over  $23^{\circ}\text{C}$  **reduces the waiting and hardening times** of the acrylic resin. Equally, **lower temperatures increase such times**.

Operation	Definition	Phase duration
Mixing	Mixing of the components	60 sec
Delivery Device Filling	The dough is transferred into the delivery device	60 sec
Waiting	The cement cannot be used	300 sec
Working	The cement can be delivered	600 sec
Hardening/Setting	The cement hardens increasing its viscosity and cannot be delivered anymore. Then the exothermic reaction takes place	360 sec

If you are using an injection system that needs to be connected to an extension tube, let a few drops of resin pour out before connecting. If the resin is injected early it will result very liquid. If the resin is injected late it will result highly viscous.

#### **THE EFFECT OF TEMPERATURE ON PREPARATION AND APPLICATION TIMING OF MENDEC SPINE**

The cement's preparation and application is strongly influenced by the temperature in the operating room. The temperature's effect on cement setting time was evaluated with a laboratory test. A graph on setting time according to temperature is provided for ease of reference (data was obtained in controlled environmental and storage conditions subjected to standard deviation). In addition to temperature and humidity, different factors can influence the cement's setting time. These include the mixing technique (speed, use of mixer), the thoroughness of mixing, the utilization of the entire liquid and solid component, the inclusion of external substances inside the cement such as blood, saline solution, etc.

#### **Caution**

- If the resin is injected early, it will result in a very liquid material; if the resin is injected late, it will result in a highly viscous material.
- The handling characteristics of bone cements are affected by operating room conditions, including the room temperature, temperature of the cement components prior to mixing, humidity, the geometry of the mixing apparatus, time spent mixing and the geometry of the delivery device. Any change in one or more of these conditions can alter the handling characteristics of bone cement.

#### **INFORMATION**

For further product information, please contact Customer Service:

Tecres S.p.A.  
Via A. Doria, 6 - 37066 Sommacampagna  
Verona Italy  
Phone: +39 0459217311  
Fax: +39 0459217330

# Mendec<sup>®</sup> Spine

High Barium Content Acrylic Resin

*for US market*

*Advancing High Technology*

TECRES S.p.A.

Via A. Doria, 6 - 37066 Sommacampagna • Verona - ITALY  
Telefono +39.045.9217311 - Fax +39.045.9217330

info@tecres.it • www.tecres.it