OP Technik
Surgical technique

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osmium is a cervical plate system for the anterior stabilization of the vertebral column. The implant immobilizes the fusion mass until the bone graft has consolidated.

The special feature of this system is the monokortical fixation with osmium screws. The osmium screw, consisting of a screw and a bolt, is expanded when the bolt is inserted.

The expansion mechanism makes superior stabilization and maximum anchorage possible, even if the bone is of reduced quality or if the same drill hole is used for a second screw placement in revision cases.

With the cancellous bone screws a bicortical fixation is achieved. The screws can be combined as per the surgeons preference.

osmium implants are made from a titanium alloy in accordance with the standard DIN EN ISO 9001.


The present surgical technique describes the working steps necessary for the application of the osmium. Also, the implants and the system instruments required for implantation are presented.

The surgical technique is not sufficient as the sole basis for the successful application of the osmium.


It is recommended to study and learn the operating technique with and from an experienced surgeon.

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The small pictures at the bottom of the page show the chronological application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-coloured.

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Indications and contraindications

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Indications

The implant system osmium is applied in cases of cervical instabilities of different origin such as conditions after anterior disc resection, fracture, tumour or pseudarthrosis resulting from previous unsuccessful operations of the cervical spine.

Contraindications

osmium must not be applied in cases of acute or chronic infection or in cases of proven or suspected intolerance to the materials used as well as in cases of poor medical or psychological general state of the patient.

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osmium screw

The osmium screw (CS 1300-14, -16, -18) is a two-component system consisting of a screw and a bolt. The implant is available in three different lengths (14, 16, 18 mm) with an outer diameter of 5 mm. The length is easily recognizable according to colour: 14 mm golden, 16 mm green and 18 mm blue. Expansion of the screw is achieved by inserting the bolt. Hyperexpansion of the screw is not possible.

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Cancellous bone screw

The cancellous bone screw (CS 1301-X) is used for the bicortical fixation. The outer diameter is 4 mm. The lengths of the screws are 12, 14, 16, 18, 20, 24 and 26 mm.

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osmium plates

The patented plates are slotted at their cranial and caudal ends. The width of the plates is 20 mm.

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Lochung	Länge
4-Loch	23, 27, 30, 33, 36 mm
6-Loch	40, 44, 48, 52, 56 mm
8-Loch	60, 66, 72 mm

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Drill guides for repositor
(two are necessary)

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Depth adjuster and
locking nut
for CS 1324-1

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Tap for cancellous bone
screws, Ø 4.0 mm,
for CS 1324

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
Surgical technique – Monocortical instrumentation

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Applying the repositor

The repositor (CS 1217-1) is positioned on the midline of the vertebral column and the drill guides (CS 1217-2) are inserted into the arms of the repositor.

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
Drilling for reduction pins

In order to fix the reduction pins (CS 1210) a drill bit of Ø 2.0 mm (CS 1215) is used.

Important: Drill under C-arm control.

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Placing of the reduction pins

After the removal of the drill guides (CS 1217-2) the reduction pins (CS 1210) are inserted into the arms of the repositor.

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Applying the reduction pins

The reduction pins are fixed by means of the screwdriver (CS 1328).

Important: Introduce the reduction pins under C-arm control. After the reduction pins have been fixed, the intersomatic space may be distracted by means of the repositor (CS 1217-1 and CS 1217-3).

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Bending of the plate

The osmium plate is contoured. Additional bending is possible if required by the patient's anatomy using the plate benders (CS 1327-1 and CS 1327-2).

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Inserting of the plate

After distraction the intersomatic space can be cleared and a bone graft or cage can be inserted. In order to insert the osmium plate the repositor must be lifted slightly.

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Compression

The achieved distraction may now be changed into compression with the lever arms for repositor (CS 1217-3).

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
Assembling the depth stop

The depth adjuster (CS 1324-2) with the locking nut (CS 1324-3) is screwed into the osmium depth stop (CS 1324-1). The drilling and tapping depth is determined by means of the x-ray images or CT/MRI images and adjusted on the depth stop. The depth stop is also used as a tissue protector. The drilling depth corresponds to the length of the screw.

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

Bore for osmium screw

The drill hole is prepared using the drill bit, Ø 2.5 mm (CS 1323) through the osmium depth stop (CS 1324). Drill under C-arm control, if necessary.

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
Tapping for osmium screw

The tap for the screw is prepared using the tap for osmium screws, Ø 5 mm (CS 1325) through the osmium depth stop.

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
Preparing for screw insertion

First remove the bolt of the osmium screw with the screwdriver (CS 1328). The attachment (CS 1310-2) is placed in the screwdriver for osmium screw (CS 1310-1) and the osmium screw is locked. The prongs of the screwdriver must be fitted to the grooves of the osmium screws.

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Implanting the osmium screws

The osmium screw is implanted.

The screwdriver is removed and the second screw on the same side is inserted.

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Implanting the osmium screws

After removal of the reposi-tor and the reduction pins, the screws are inserted at the opposite side in the same way. Using the screwdriver (CS 1310-1) or the tightener for osmium screws (CS 1322) the screws can be further tightened before being expanded, in order to achieve an even support of the wings.

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
Fixing the osmium screws

The screwdriver for osmium screw (CS 1310-1) without the attachment (CS 1310-2) is placed onto the osmium screw and the bolt is inserted into the shaft of the instrument. The bolt is locked using the screwdriver, hex 2.5 mm (CS 1328) and the osmium screw is expanded.

Important: The screwdriver for osmium screw must remain engaged on the grooves of the osmium screw, in order to prevent that the screw co-rotates.

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Fixation of the bone graft

If necessary the bone graft can be fixed using a cancellous bone screw in the whole in the center of the plate.

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Surgical technique – Bicortical instrumentation

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Mounting the repositor, insertion of the plate, and compression

Use the same technique as with monocortical screws.

Preparing the screw implantation

For the cancellous bone screws the same drill bit (CS 1322) and the depth stop (CS 1324) may be used. Determining the drill depth by using the depth stop is possible as well. The drilling depth corresponds to the screw length.

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Tapping for cancellous screws

Use the same technique as with the monocortical screws. However, use the tap for cancellous bone screws, Ø 4.0 mm (CS 1326).

Important: It is optional to use the tap (CS 1326) before applying the screws. Use a C-arm.

Implanting the cancellous bone screws

The appropriate length of the cancellous bone screws may be determined using the depth gauge (UT 1068-23). Inserting the cancellous bone screws by means of the screwdriver (CS 1328) under C-arm control.

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







Components

Implants	Product number
osmium screw incl. bolt, golden, Ø 5 mm, titanium, length 14 mm	CS 1300-14T
osmium screw incl. bolt, green, Ø 5 mm, titanium, length 16 mm	CS 1300-16T
osmium screw incl. bolt, blue, Ø 5 mm, titanium, length 18 mm	CS 1300-18T
Cancellous bone screw, Ø 4 mm, titanium, length 12 mm	CS 1301-12T
Cancellous bone screw, Ø 4 mm, titanium, length 14 mm	CS 1301-14T
Cancellous bone screw, Ø 4 mm, titanium, length 16 mm	CS 1301-16T
Cancellous bone screw, Ø 4 mm, titanium, length 18 mm	CS 1301-18T
Cancellous bone screw, Ø 4 mm, titanium, length 20 mm	CS 1301-20T
Cancellous bone screw, Ø 4 mm, titanium, length 22 mm	CS 1301-22T
Cancellous bone screw, Ø 4 mm, titanium, length 24 mm	CS 1301-24T
Cancellous bone screw, Ø 4 mm, titanium, length 26 mm	CS 1301-26T
osmium plate, 4 holes, titanium, width 20 mm, length 23 mm	CS 1306-23T
osmium plate, 4 holes, titanium, width 20 mm, length 27 mm	CS 1306-27T
osmium plate, 4 holes, titanium, width 20 mm, length 30 mm	CS 1306-30T
osmium plate, 4 holes, titanium, width 20 mm, length 33 mm	CS 1306-33T
osmium plate, 4 holes, titanium, width 20 mm, length 36 mm	CS 1306-36T
osmium plate, 6 holes, titanium, width 20 mm, length 40 mm	CS 1307-40T
osmium plate, 6 holes, titanium, width 20 mm, length 44 mm	CS 1307-44T
osmium plate, 6 holes, titanium, width 20 mm, length 48 mm	CS 1307-48T
osmium plate, 6 holes, titanium, width 20 mm, length 52 mm	CS 1307-52T
osmium plate, 6 holes, titanium, width 20 mm, length 56 mm	CS 1307-56T
osmium plate, 8 holes, titanium, width 20 mm, length 60 mm	CS 1308-60T
osmium plate, 8 holes, titanium, width 20 mm, length 66 mm	CS 1308-66T
osmium plate, 8 holes, titanium, width 20 mm, length 72 mm	CS 1308-72T
Instruments	Product number
Reduction pin	CS 1210
Drill bit for reduction pin, Ø 2.0 mm	CS 1215
Repositor	CS 1217-1
Drill guides for repositor	CS 1217-2
Lever arm for repositor	CS 1217-3
Screw holding forceps	CS 1219
Screwdriver for osmium screw	CS 1310-1
Attachment for screwdriver	CS 1310-2
Plate holding forceps	CS 1319
Bone graft holder	CS 1320
Bone graft impactor	CS 1321
osmium screw tightener	CS 1322
osmium drill bit, Ø 2.5 mm, for CS 1324	CS 1323
osmium depth stop	CS 1324-1
Depth adjuster for CS 1324-1	CS 1324-2
Locking nut for CS 1324-1	CS 1324-3
Tap for osmium screws, Ø 5 mm, for CS 1324	CS 1325
Tap for cancellous bone screws, Ø 4.0 mm, for CS 1324	CS 1326
Plate bender	CS 1327-1
Plate hold and bend instrument	CS 1327-2
Screwdriver, hex 2.5 mm	CS 1328
Measuring caliper	CS 5788
Depth gauge, length 230 mm	UT 1068-23

Ulrich GmbH & Co. KG

Buchbrunnenweg 12

89081 Ulm

Germany

Telefon/Phone +49 (0)731 9654-110/225

Telefax/Fax +49 (0)731 9654-2702

e-mail spine@ulrichmedical.com

internet www.ulrichmedical.com