

# Supplement für OP-Technik Supplement for Surgical Technique

## uCentum™

comprehensive posterior system

**Perkutan**  
**Percutaneous**

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+ Percutaneous distraction  
and compression

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The small pictures at the bottom of the page show the step-by-step application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-colored.

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Retractor for percutaneous  
reduction, left

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Retractor for percutaneous  
reduction, right

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Retractor for percutaneous  
reduction, variable

**CS 3843-04**

Repositionsarm,  
perkutan, links

Reduction arm,  
percutaneous, left

**CS 3843-05**

Repositionsarm,  
perkutan, rechts

Reduction arm,  
percutaneous, right

**CS 3843-06**

Repositionsarm,  
perkutan, variabel

Reduction arm,  
percutaneous, variable

**CS 3843-07**

Repositionsrad

Reduction wheel

**CS 3843-08**

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Nut for retractor,  
percutaneous

**CS 3843-20**

Distraktionszange,  
perkutan



Distraction forceps,  
percutaneous

**CS 3843-21**


Kompressionszange,  
perkutan

Compression forceps,  
percutaneous


✓ Erlaubte Distraktion und Kompression ✗ Verbotene Distraktion und Kompression		1 Zange	2 Zangen	1 Sperrer (V)	2 Sperrer (V+L oder R)
✓ Allowed distraction and compression ✗ Prohibited distraction and compression		1 forceps	2 forceps	1 retractor (V)	2 retractors (V+L or R)
					
CS 3802-xx-xxx	Schraube, polyaxial Screw, polyaxial	✓	✓	✓	✓
CS 3803-xx-xxx	Schraube, monoaxial Screw, monoaxial	✗	✗	✗	✗
CS 3804-xx-xxx	Schraube, polyaxial, vorklemmbar Screw, polyaxial, pre-fixation	✓*	✓*	✓*	✓*

[1] Winkelvariable Distraktion und Kompression  
[1] Variable-angle distraction and compression

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[2] Winkelvariable Distraktion  
[2] Variable-angle distraction


[3] Winkelvariable Kompression  
[3] Variable-angle compression

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With the components of the percutaneous distraction/compression instrument set, anatomic profiles can be percutaneously reconstructed and secured through distraction or compression via the positioned rod. To be able to perform this application in this way, the rod must be finally secured in one screw head. Then the distraction or compression can be performed against the guide sleeve of the secured screw.

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
As a general rule, a differentiation is made between variable-angle [1 to 3, open application here] and stable-angle [4 to 6, open application here] distraction/compression, whereby polyaxial screws are used for the variable-angle [1] (lordotic or kyphotic) application and monoaxial or pre-fixed screws are used for the stable-angle [4] (parallel) distraction/compression.


✓ Erlaubte Distraktion und Kompression ✗ Verbotene Distraktion und Kompression ✓ Allowed distraction and compression ✗ Prohibited distraction and compression		1 Zange	2 Zangen	1 Sperrer (V)	2 Sperrer (V+L oder R)
		1 forceps	2 forceps	1 retractor (V)	2 retractors (V+L or R)
					
CS 3802-xx-xxx	Schraube, polyaxial Screw, polyaxial	✗	✗	✗	✗
CS 3803-xx-xxx	Schraube, monoaxial Screw, monoaxial	✓	✓	✓	✓
CS 3804-xx-xxx	Schraube, polyaxial, vorklemmbar Screw, polyaxial, pre-fixation	✓*	✗*	✓*	✗*

[4] Winkelstabile Distraktion und Kompression  
 [4] Stable-angle distraction and compression

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There are four techniques available:

- Distraction/compression with 1 forceps
- Distraction/compression with 2 forceps
- Distraction/compression with 1 retractor
- Distraction/compression with 2 retractors

**Caution:** Only straight or slightly curved rods may be used.

**Important:** Perform the techniques described below under continuous C-arm control.

**Important:** In all of the techniques described below, the screw head can become jammed on the rod, particularly when the two-part retractor is used and in the case of long reduction channels. Through manual application of a small amount of force in the opposite direction on the upper end of the sleeve, this tension can be loosened and mobility along the rod can be restored.

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**Important:** The locking screw in the opposite screw head is already inserted but has not yet undergone final tightening. When using two pairs of forceps or the 2-part retractor the second locking screw must not be inserted during reduction.

The final fixation of the second locking screw takes place following reduction and after disassembling the selected reduction instruments with the help of the mounted counterpart, percutaneous (CS 3849 + CS 7144) and the corresponding T-handle with torque limitation (corresponding to the selected screw).

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### ■ 1. Distraction / compression with one pair of forceps

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**Important:** In both techniques described below, particularly in the case of longer reduction channels, the screw head may become jammed on the rod. Through manual application of a small amount of force in the opposite direction on the upper end of the sleeve, this tension can be loosened and mobility along the rod can be restored.

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
#### Distraction with one pair of forceps

For distraction of both vertebral bodies, the distraction forceps, percutaneous (CS 3843-20) is positioned between the guide sleeves as close to the level of the skin as possible and spread apart until the semicircular adapters are in direct contact with the guide sleeves. The two vertebral bodies are now pushed apart over the screw heads by spreading the distraction forceps [1] further apart.

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#### Compression with one pair of forceps

For compression of both vertebral bodies, the compression forceps, percutaneous (CS 3843-21) is positioned from the outside on the guide sleeves as close to the level of the skin as possible and pressed together until the semicircular adapters are in direct contact with the guide sleeves. The two vertebral bodies are now pressed together over the screw heads by pressing the compression forceps [2].

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
**Caution:** When using the two pairs of forceps, the locking screw should not yet be screwed into the opposite screw head in order to avoid tensioning of the overall construction. As a result, the pre-fixation feature of the screws, pre-fixation (CS 3804-xxx-xxx) is not applicable in case of two forceps are used.

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
**Distraction with two pairs of forceps**

First the distraction forceps, percutaneous (CS 3843-20) are placed as far down as possible at the level of the skin and spread apart until the semicircular adapters are in direct contact with the guide sleeves. Subsequently the compression forceps, percutaneous (CS 3843-21) are placed as far up as possible on the guide sleeves. The distraction of both vertebral bodies then takes place by pressing both guide sleeves together further with the compression forceps. The distraction forceps serve as a counterpart [3] in this case.

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
**Compression with two pairs of forceps**

First the compression forceps, percutaneous (CS 3843-21) are placed as far down as possible at the level of the skin and pressed together until the semicircular adapters are in direct contact with the guide sleeves. Subsequently the distraction forceps, percutaneous (CS 3843-20) are then placed as far up as possible between the guide sleeves and spread apart. The compression of both vertebral bodies then takes place by spreading both guide sleeves apart further with the distraction forceps. The compression forceps serve as a counterbearing [4] in this case.

### ■ 3. Montage des Repositionsinstruments

### ■ 3. Assembly of the Reduction Instrument

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The assembly is performed in the same way for all three instruments (R, L, V). It is shown here as an example.

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
**Important:** The button on the reduction arm must be snapped into place vertically for assembly [5].

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


**Important:** It should be ensured that the reduction arm and retractor with the same letters and the same alignment are pre-assembled together [6 and 8].

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The reduction wheel (CS 3843-07) is initially inserted in the reduction arm (CS 3843-04/-05/-06) [6]. The preassembled reduction arm is pushed onto the retractor (CS 3843-01/-02/-03) with the head positioned vertically [7 and 8].

■ 4. Distraction / Kompression mit einem Sperrer  
■ 4. Distraction / compression with one retractor

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### Assembly of the Reduction Instrument on the Guide Sleeves

The assembled variable reduction arm **V** is guided over both guide sleeves in place with its head positioned vertically and pushed down to the skin [9].

**Important:** The variable reduction arm must be placed as low down as possible at the level of the skin.

The distraction and compression may now be performed either by hand or with the aid of the screwdriver for rod inserter, SR 20 (CS 3831). Compression or distraction (according to the direction of the arrow) can be set by pulling and then rotating the button on the reduction arm [10].

**Important:** Particularly in the case of longer reduction channels, the screw head may become jammed on the rod. Through manual application of a small amount of force in the opposite direction on the upper end of the sleeve, this tension can be loosened and mobility along the rod can be restored.

Following the reduction, the position of the rod must be right-angled to the screw head under C-arm control or corrected, if applicable [11]. Prior to final fixation with the corresponding torque, the variable retractor must be removed [12].

## ■ 5. Distraction / Kompression mit zweiteiligem Sperrer ■ 5. Distraction / compression with a two-part retractor

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
For reduction with a two-part retractor, two mounted reduction arms are needed for each side. For each arm on top with fixed orientation (L or R), an arm on the bottom with variable orientation is needed (V). Prior to assembly, the surgeon determines the caudal/cranial position of the toothed rods in the surgical field. If required, the two-part retractor can also be mounted on both sides.

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**Caution:** When using the two-part retractor, the locking screw should not yet be screwed into the opposite screw head in order to avoid tensioning of the overall construction. As a result, the pre-fixation feature of the screws, pre-fixation (CS 3804-xxx-xxx) is not applicable in case of two forceps are used.

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
### Assembly of the Reduction Instrument on the Guide Sleeves

The assembled variable reduction arm **V** is guided over both guide sleeves in place with its head positioned vertically and pushed down to the skin [13].

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Then the reduction instrument (**L** or **R**) in the opposite direction is placed on the guide sleeves [14].

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




To secure the distal reduction instrument, the nut for retractor, percutaneous (CS 3843-08) must be screwed in as far as it will go [15 and 16].

**Important:** It should be ensured that both reduction wheels (CS 3843-07) can be operated from above.

**Important:** The procedure must be performed under continuous C-arm control.

Compression or distraction (according to the direction of the arrow) can be set by pulling and then rotating the button on the reduction arm [17].

The distraction and compression may now be performed either by hand or with the aid of the screwdriver for rod inserter, SR 20 (CS 3831) [18]. Reduction to the desired position is performed by alternately rotating both reduction wheels.

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**Important:** The rod must definitively be right-angled to the screw head before the locking screw can be introduced [19]. If necessary, the position of the screw head can be corrected by the reduction wheels. If the rod is not completely in the screw head, the rod must be completely inserted in the screw head by using the rod inserter with spring, percutaneous (CS 3829) or, if pre-fixation screws are used, the rod inserter CS 3829-01.

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
After reduction and position control of the rod in the head have been performed, the second locking screw can also be introduced. Before final fixation, the upper reduction arm must be unlocked in the neutral position by loosening the nuts (CS 3843-08) as well as by pulling and rotating the button on the reduction arm and the reduction instrument must be removed from the sleeves.

**Caution:** During loosening, there may be a slight loss of correction. As a result, however, the screw head aligns itself vertically to the rod which is absolutely necessary for secure insertion of the locking screw. By rotating the lower reduction wheel, the loss of correction can be compensated.

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The second locking screw must now be initially tightened by hand. Then the lower variable retractor is detached by unlocking the button on the reduction instrument in the neutral position by pulling and rotating and the instrument is removed from the sleeves. Then the locking screw with the corresponding torque undergoes final fixation [20].

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**Ulrich GmbH & Co. KG**