


# ALIF Footprint. Less Invasive. Fusion Without Compromise.™

**InterFuse® Interbody Fusion System Featuring the InterFuse S™ Procedure**

**Achieve an ALIF-Sized Footprint and Utilize a Familiar Posterior Approach**



- **Insert** module down rail
- **Move** modules laterally with Positioning Lever
- **Remove** rail on previous module
- **Repeat**

(1-2) After gaining access to the spine via a less invasive unilateral posterior approach, the initial module is placed into the prepared disc space. (3-5) Additional modules are inserted into position via InterFuse's unique rail and slot technology. (6) InterFuse's modular, intra-operative assembly™ provides a large footprint that can be customized (3-6 modules) for each patient.