

# JAZZ™ PASSER

JAZZ™ Passer is an additional solution dedicated to spine surgeons bringing them a new comfort in using Band Connectors.

JAZZ™ Passer is a new system:

- 3 types of Passers,
- 1 Band-sublaminar, subtransverse, interspinous.

## Passers

- **Straight**



- **Round**



- **Round-Angled**



## JAZZ™ Passer Band

**Reduced thickness** with a distal hole which is grasped with the JAZZ™ Passer hook.

## Universal

JAZZ™ Passer can be used with all the JAZZ™ Platform Band Connectors.



JAZZ™ Lock



JAZZ™ Standalone



JAZZ™ Band



JAZZ™ Claw

## Simple

A simple rotation of the Passer handle passes the Band around the anatomical structure.

No additional instrument tray required.

The 3 types of Passers integrated into existing JAZZ™ Band trays.



# REFERENCES



Straight Passer

Reference	Designation
550341	Straight Passer - 15mm
550342	Straight Passer - 20mm
550343	Straight Passer - 25mm



Round Passer

Reference	Designation
550344	Round Passer



Round Passer-Angled

Reference	Designation
550345	Round Passer-Angled



JAZZ™ Passer Band

Reference	Designation
150157	JAZZ Passer Band

JAZZ™ Passer Band is used with the JAZZ™ Platform Connectors