



PRECISION SPINE

SURELOK™ MIS 3L

PERCUTANEOUS SCREW SYSTEM

SIMPLE & EFFICIENT

The SureLOK™ MIS 3L is a low profile, percutaneous screw system. The system features an extended tab tulip with 20mm of controlled rod reduction. The screw has a triple lead, proximally tapered thread that greatly enhances delivery. This reduces fatigue and minimizes OR time without compromising screw purchase. In addition, the open tulip design, 30° of rod angulation and rod tip eases rod insertion. The system also offers multiple sleeve options to allow for aggressive tulip manipulation, rod reduction and compression and distraction techniques.



SURE-LOK™ MIS 3L PERCUTANEOUS PEDICLE SCREW

150mm Integrated Extension Tab features a low profile outer diameter (13.7mm) allowing for a smaller incision to minimize muscle disruption

Open Tab Design provides maximum visualization and eases rod placement

Proximal Window allows for rings and sleeves to securely lock into place

Polyaxial Screws offer 30° of angulation in all planes for intraoperative flexibility

Integrated Reduction Threads provide 20mm of controlled rod reduction


Aggressive Self Drilling, Self Tapping Screw Tip optimizes screw engagement

Break-Off Tab removal is completed with one easy step


Triple Lead Thread allows for efficient screw delivery

Proximal Tapered Thread design increases bone screw interface, enhancing pull-out strength

Curved & Straight Rod's bullet-tipped end and keyed hex end provide a rigid connection while the 15° rod inserter angle eases rod insertion and release



Open Tulip Design, 30° of rod angulation and rod tip eases rod insertion



Streamlined Insertion Instrumentation provides OR efficiency while minimizing soft tissue disruption:

- **Dilator/Awl/Tap** is completed with one instrument



- **Modular Rings & Sleeves** provide maximum rigidity for rod manipulation



- **External Rod Reducer** allows for controlled rod reduction



- **Compressor/Distractor** is powerful, precise and effective



MIS 3L Rod Sizing – 5.5mm ø

Ti	Ti
Straight	Curved
35 - 80mm	35-80mm (5mm)
90mm	90-150mm (10mm)
100mm	
110mm	

MIS 3L Sizing

5.5mm Diameter	6.5mm Diameter	7.5mm Diameter	8.5mm Diameter
35mm	35mm	35mm	35mm
40mm	40mm	40mm	40mm
45mm	45mm	45mm	45mm
50mm	50mm	50mm	50mm
55mm	55mm	55mm	55mm

Precision Spine, Inc.

2050 Executive Drive, Pearl, MS 39208
 Customer Service: 1.888.241.4773
 Phone: 601.420.4244
 Toll Free: 877.780.4370
 Fax: 601.420.5501
www.precisionspineinc.com



CE0086